

Technical Data

Product: THERMOVAC TM 101
Part No.: 230081V01

| | |
|----------------------------------|---|
| Vacuum connection: | G1/4 female thread / DN 16 ISO-KF adapter |
| Measurement range: | 1200 to 5×10^{-4} mbar |
| Measuring units: | mbar, Torr, Pa |
| Display: | LCD 12 mm |
| Measuring cycle: | 1.6 s |
| Measuring accuracy | |
| at 1200 - 10 mbar: | 0.3% of full scale |
| at $10 - 2 \times 10^{-3}$ mbar: | 10 % of measured value |
| $< 2 \times 10^{-3}$ mbar: | faktor 2 |
| Admissible ambient temperature: | +5 to +50°C |
| Power supply: | 9 V block battery / Plug-in mains adapter |
| Max. operating time: | 75 h |
| Protection class: | IP 40 |
| Dimensions: | see dimension sheet |
| Weight: | 0.23 kg |

Technical data are subject to change